



vh 2100

Humidifier

- ⚠ Designed to deliver a comfortable level of humidity.

The VH-2100 Humidifier is part of a complete humidification system designed to deliver a comfortable level of humidity. Adjustable between nine temperature and humidification levels.

This arrangement simplifies setup, minimizing the connections required.

Product specifications

	VH-2100-115	VH-2100-230
Dimensions	W135 x D170 x H156 mm	
Weight	1.5 Kg (w/o Chamber fitted)	
ELECTRICAL POWER		
Supply Voltage	115 ± 12V	230 ± 23V
Supply Current	1.3 A Max.	0.7 A Max.
Supply Frequency	50/60 Hz	
Heater plate capacity	150W max	
SECURITY PROTECTION		
Overcurrent protection	115 V/3 A (fast burn type)	230 V/1.5 A (fast burn type)
Overheat protection	115±3° C, thermostat.	
TEMPERATURE SETTING		
	Knob: Scale 1 to scale 9	

55+ years of experience committed to innovation and patient care.

www.tecmeglobal.com / info@tecmeglobal.com

Córdoba - Buenos Aires - Atlanta - São Paulo